



EYE/FACE WASH & BODY SPRAY

LSE 1B

G1 EYE/FACE WASH PUSH LEVER OPERATED – WALL MOUNTED

RECEPTOR BOWL

High Visibility Yellow extra-large Cicolac receptor Bowl resistant to damage from alkalies, salt solutions, oils, most acids.

ATOMIZER HEADS

Twin-cushion-flow Cicolac Atomizer Heads in Safety Green and Yellow with flip-flop dust covers that protect heads and automatically releases with water pressure.

ACTIVATION

Yellow Cicolac Push plate/Foot pedal operated Stay-open Brass Ball valve with Stainless Steel Ball and Spindle.

STREAM CONTROL

Steady water flow under varying pressure assured with Automatic pressure compensators for 30-90 PSI.

SUPPLY

1/2" BSP

Recommended pressure: 30 PSI/2 kg./ cm²

DRAIN

1 1/4" BSP

NET SHIPPING WEIGHT

3 kg

